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Instructions for Use

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MADISON Total knee
prosthesis

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2010

You have chosen an implant intended for knee arthroplasty manufactured by IMPLANET. As such it satisfies the requirements relating to design, manufacture and control described in standard ISO 13485.

- Prior to use it is essential to read through the information provided in these instructions as well as the data shown on the labeling.
- If applicable, the accompanying documents such as sales documentation, operating technique, tracings, technical files and kit follow-up plans are available from our customer service and/or on the IMPLANET internet site.
- This implant must be used exclusively by trained and qualified persons who are familiar with these instructions for use.

1. Performance

- This implant forms part of a range designed and produced in order to restore the functionality of the knee joint within the context of first intention treatment or in the event of a revision.
- The longevity and the quality of the arthroplasty depend on its environment, the positioning conditions and its conditions of use.
- In order to obtain optimum performances and longevity of the prosthesis, it is essential to comply with the combinations recommended in our documentary system.
- Under no circumstances must implants from the MADISON range be mixed with implants from another range.

2. Materials

- IMPLANET uses exclusively, unless specified otherwise, the following standardized materials specifically intended for surgical implants:

Femoral components and tibial base plates:

- Chrome/cobalt alloy
- Titanium alloy
- Pure titanium and hydroxyapatite (coated implants only)

Patellas:

- Very high molecular weight polyethylene
- Titanium alloy

Inserts:

- Very high molecular weight polyethylene
- Uncoated femoral components and tibial base plates are fixed with acrylic cement.
- The exact nature of the material of each component is mentioned on the label.
- As far as assembly is concerned IMPLANET respects the compatibility of the materials defined by the regulations in force.

3. Indications

This implant is intended to replace the knee joint and the main indications are as follows:

- Degenerative osteoarthritis or post traumatic arthrosis, in patients whose age, weight, and level of activity are compatible with a satisfactory long term result
- Rheumatoid polyarthritis
- Revision of unicompartmental arthroplasty procedure, if bone mass and ligament stability allow
- Failed osteotomy

4. Contraindications

It is usual to consider some of these problems as contraindications to positioning a total knee prosthesis:

Absolute contraindications:

- Inadequate or poor quality bone support,
- lateral ligaments damaged, malfunction of extensor apparatus,
- local infection, some time ago or earlier, sepsis of the knee,
- serious vascular disease,
- recurvatum deformity originating in muscle deficiency,
- flexum deformity which cannot be corrected,
- successful arthrodesis of the knee,
- parameters which are incompatible with satisfactory long term results: age, weight, level of physical activity, etc.,
- neurological condition of the patient incompatible with the post-operative constraints associated with this type of intervention.

Relative contraindications:

- medical contraindication to anesthesia,
- known allergy to one of the implant components,
- history of osteomyelitis around the knee,
- neuropathic joint,
- muscle or tissue cover of the operational site inadequate,
- obesity.

However, the surgeon is responsible for evaluating the risk/benefit ratio when making a decision about whether or not to use a prosthesis within the context of these etiologies.

N.B.: when positioning an implant it is considered to be the best solution for the patient and the latter presents one or several contraindications. The surgeon is responsible for informing this patient of the possible risks regarding the success of the intervention and for giving him or her any useful advice relating to the measures which may reduce the effects of these contraindications.

5. Storage conditions/ Expiry date

- The storage and handling conditions must be able to guarantee the integrity of the implant and its packaging.
- Damaged implants must not be used.
- Before using this implant, it is essential to make sure that the packaging is undamaged and to check the expiry date indicated on the label, which provide a guarantee that the implant is in a sterile condition.

6. Sterilization

- This implant was sterilized in accordance with current regulations by gamma irradiation ("R" labeling).
- Re-sterilization of this implant is categorically prohibited and would free IMPLANET from any liability.
- If a single use product is re-used, the performance, cleaning and sterility of the device can no longer be guaranteed. This can lead in particular to failure of the intervention or to the risk of infection which may result in the patient's death.

7. Recommendations for use

- It is essential to comply with conditions of asepsis when opening the protective packaging and removing the implant.

- Certain parts of the implant are specifically protected or are delivered with a single use device intended to facilitate their use. Take care when removing this protection or these devices once they are no longer useful.
- All implants must be used in their original form unless specifically mentioned otherwise. Any required changes to the implant are the surgeon's exclusive responsibility. Implants must be handled with meticulous care; cutting, twisting or scratching their surface may significantly reduce their lifetime.
- There is a right and left version of the femoral component. It is essential to comply with the indication for the appropriate side clearly indicated on the label.
- The size of the permanent implants can be determined via pre-operative planning and use of trial implants.
- Only correct use of the ancillary equipment supplied by IMPLANET and specific to the implant guarantees its satisfactory positioning.
- In the event of explantation, decontamination and sterilization conditions must be respected.
- We categorically advise against combining our products with elements from other sources (manufacturer).
- The waste resulting from the operation (packaging, explants, etc.) must be dealt with in the same way as the health facility deals with any other medical waste.
- The cementing technique requires special precautions (preparation of the bone cavity, cementing technique, cooperation with the anesthetist, etc.) which are listed in the special instruction leaflet for the cement.
- MADISON devices have not been the subject of a safety and compatibility assessment in a magnetic resonance environment. Heating and migration of MADISON devices in a magnetic resonance environment have not been tested.

7.1 Recommendations for general use relating to technique

- Surgical actions aimed at preparing the bone bed must be adapted to suit the dimensions and mechanical qualities of the bone.
- The stability and positioning of the modular assemblies must always be checked.
- The choice of approach to the knee and its closure must take into account the quality of the architectural and mechanical reconstruction of the knee and must respect muscle balance.
- In all cases the quality of the implant (stability, positioning, cleanliness of surfaces, search for cam effect, reduction test, etc.) must be checked once positioning is complete.

7.2 Specific recommendations

- The 3 main elements constituting the total prosthesis of the knee (femoral component, insert, tibial base plate) are each available in 8 sizes, from 1 to 8. Dimensional differences are possible between the femoral condyles, on the one hand, and the tibial base plate, on the other hand, in the same patient. Consequently the surgeon may prefer to choose a femoral component and a tibial base plate of different sizes.

Generally speaking, the MADISON prosthesis tolerates a difference of 2 sizes when the tibial component is smaller, and with a mobile plate, and a difference of more or less one size in all other configurations. Implanet recommends not using components with differences greater than these.

There are two possible cases:

- a prosthesis with a mobile plate. For differences of one size, the insert can be the size of either the tibial component or the femoral component. For differences of 2 sizes, the insert must be of the intermediate size. Example: if a size 3 tibial plate and a size 5 femur are used, the insert must be size 4.
- a prosthesis with a fixed plate; the tibial base plate and the insert must be of the same size.

Mobile plate			
	Difference of 1 size between femoral component and tibial base plate		Difference of 2 sizes between femoral component and tibial base plate
	1st case: femoral component and insert are the same size	2nd case: insert and tibial base plate are the same size	
Femoral component	0	Size +1 or -1	Size +1
Insert		0	0
Tibial base plate	Size +1 or -1		Size -1

Fixed plate		
	Difference of 1 size between femoral component and tibial base plate	Not applicable
Femoral component	Size + 1 or -1	
Insert	0	
Tibial base plate		

- The choice of either cementing one of the two metal components (femur or tibial base plate) or not doing so is independent of the choice of the other component. We do not have any data concerning combinations of a cemented femur with an uncemented socket or an uncemented femur with a cemented socket: use of these configurations is contraindicated.
- If there are any doubts about the stability of the definitive tibial base plate, extension stems are provided and they are positioned by screwing under the tibial keel, after removing the pre-assembled metal cap.
- The designations ML and AP on the labels indicate the Medio-Lateral and Antero-Posterior dimensions of the implants.

8. Adverse effects

8.1 Adverse side effects relating to the patient:

- deep vein thrombosis and pulmonary embolism,
- hypersensitivity reaction to one of the implant components.

8.2 Adverse side effects associated with arthroplasty:

- bone fractures (including fracture of the patella),
- reactions in bones in contact with the implant, osteolysis,
- articular or peri-articular calcifications,
- joint failure: failure of the extensor apparatus, instability and/or limitation of joint mobility associated with poor positioning of the implant, flexum or recurvatum, faulty alignment of the knee, varus or valgus deformities, excessive rotation,

- stiffness, luxation or subluxation,
- unequal limb length,
- scars or delayed healing,
- hematoma, vascular lesions and/or thromboembolic accident,
- peripheral neuropathies,
- infection,
- pain,
- loosening which may cause migration, deformity, cracks or fractures of a component,
- failure of hydroxyapatite / titanium spray adhesion,
- insert wearing at the joint surfaces where friction occurs,
- local or remote tissue reaction caused by debris from wear inherent in the functioning of a knee prosthesis (this debris may spread inside the body and start up a tissue reaction),
- allergic sensitivity reaction to the materials constituting the implants,
- using cement has its own complications.

9. Information to be given to patients by the surgeon

- The surgeon must advise patients with implants to comply with the rules relating to a hygienic lifestyle and warn them of any incorrect or harmful stress applied to the joint with the prosthesis. It is particularly advisable to inform patients that their weight and level of activity may influence the longevity of the prosthesis.

Factors likely to compromise the success of the implant operation are, specifically:

- non-prevention of local and general infections,
 - drug addiction and/or habit of drug and medication abuse,
 - excess weight,
 - intense physical activities or repeated movements which expose the prosthesis to excessive loads (hard physical work, running, etc.),
 - risks of conflict with other prostheses,
 - patient's intellectual incapacity to understand and follow the doctor's instructions.
- At the time of any treatment or investigation in the vicinity of the implant (injections of medication into the buttocks, or scan, or MRI, etc.), the patient must indicate that he has a prosthesis.
 - The surgeon must inform the patient that as time goes by there is a possibility of experiencing wear of certain prosthetic components and possible changes in the bone in contact with the implant. He or she must be informed that regular

monitoring of the prosthesis by an orthopedic surgeon is essential and he or she must be advised to arrange another consultation for any symptoms which seem abnormal.

10. Liability

IMPLANET cannot be held responsible for any consecutive or accidental loss relating to the operation, damages or costs, directly or indirectly resulting from using the product. IMPLANET does not assume and does not authorise any third party to take on in its name any additional liability nor any liability associated with its products.

11. Information on quality

If this implant fails to meet with your entire satisfaction, please contact the IMPLANET quality department and notify us of the nature of the faults noted with the reference and the lot number for the implant.

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